Model **23650-Series** POLYPROPYLENE

LABORATORY **DRUM PUMPS**

Tube:

Polypropylene

Impeller:

Axial Flow, Polypropylene

Shaft:

Hastelloy C Lip Seal of Viton*

Seal: Port:

3/4" Dia., Slip-on Hose Discharge

Weight:

1- 3/4 lb. (0.8 Kg) average



Specifically designed to be used in laboratory size containers where fluids are to be transferred to smaller containers for handling.



Explosion hazard. Do not use where flammable vapors are present. Do not pump solvents, thinners or gasoline. To do so can cause an explosion which can result in injury or death.

OPERATING INSTRUCTIONS

- General When operating the JABSCO drum pump, protective clothing, gloves, and safety glasses should be worn. All standard safety practices should be followed.
- 2. Installation Pump may be mounted through bung hole in drums, on side of vats and tanks, and through top of carboys. Barrel Adaptors are available to hold pump in stationary or rigid position. When inserting drum pump into barrel, vat or drum slowly lower pump into container to avoid spills. Correct pump length should be selected for container with which pump is to be used.
- Drive An electric motor (23640-series) and a rechargeable motor (23780-series) are available to drive the pump. Refer to Motor Data Sheet for details (form number 43000-0690 and 43000-0689).
- **Temperature Range** Maximum pumpage temperatures are 32°F 130° F (0° 55° C).
- Performance Limits Maximum specific gravity of 1.8 or viscosity of 2500 ssu (540 cp) are the recommended limits for normal use. Contact factory for further information.
- Spare Parts Use JABSCO service parts for replacement and
- 7. Running Dry Unit has been designed for limited run dry capacity. Do not run dry for more than 30 minutes. Lack of liquid will eventually damage bearing and other pump parts.
- * Viton is a registered trademark of DuPont Dow Elastomers L.LC.

Jabsco



U.S.A. ITT Jabsco 1485 Dale Way, P.O. Box 2158 Costa Mesa, CA 92628-2158 Tel: (714) 545-8251 Fax (714) 957-0609



UNITED KINGDOM ITT Jabsco Bingley Road, Hoddesdon Hertfordshire EN11 OBU Tel: +44 (0) 1992 450145 Fax: +44 (0) 1992 467132

CANADA ITT Fluid Products 55 Royal Road Guelph, Ontario N1H 1T1 Tel: (519) 821-1900 Fax: (519) 821-2569

JAPAN NHK Jabsco Company Ltd. 3-21-10, Shin-Yokohama Kohoku-Ku, Yokohama, 222 Tel: 045-4758906 Fax: 045-475-8908

GERMANY Jabsco GmbH Oststrasse 28 22844 Norderstedt Tel: +49-40-53 53 73 -0


Rev. 9/1998

Model 23650-Series

VARIATIONS AVAILABLE

MODEL NO.	DESCRIPTION
23650-0500	20" (500 mm) POLYPROPYLENE TUBE
23650-0700	28" (700 mm) POLYPROPYLENE TUBE
23650-1000	39 1/2" (1000 mm) POLYPROPYLENE TUBE

PERFORMANCE TABLE



FLOW - LITERS PER MINUTE

THE PRODUCTS DESCRIBED HEREIN ARE SUBJECT TO

IS AVAILABLE FOR YOUR INSPECTION UPON REQUEST.

THE JABSCO ONE YEAR LIMITED WARRANTY, WHICH

Form: 43000-0687

© Copyright 1998, ITT Industries

Printed in U.S.A.

All Rights Reserved

Fax: +49-40-53 53 73 -11